

Medical CarePulp® : Versatile Solutions

BULBOUS URINAL

Registered medical product with 



875ML

Dimensions piece 260mm x 130mm x 100mm (L x W x H)

Weight / 30 g

piece Carton

Weight of the Carton : 30 g

Carton dimensions: 260 mm x 100 mm x 130 mm (L x W x H)

Number of pieces per carton : 100

Number of carton per palette : 5

Number of pieces per pallet : 500

Pallet dimensions: 1200mm x 1000mm x 1028mm (L x W x H)

No printing

Available colors

Grey

UDI-DI 5060802530001

Product should be kept in a dry environment.
If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).
The characteristics are given as an indication. Non contractual document. Changes excepted. Weights & measures are of average data!
You can find more information about this product on the following page

Version 13.06.2023



Omni-Pac Group
Second Avenue,
Flixborough Industrial Estate
Scunthorpe
North Lincolnshire DN15 8SD
UNITED KINGDOM
Tél: +44 (0)1724 847 207
www.omnipacgroup.co.uk



Supplement to the product data sheet

Medical products

Product description

Composition (materials used):

Moulded fibre products made of recycled paper and fresh fibre or only recycled paper with the addition of process-related and functional additives.

Intended use:

As medical devices for use in the medical field, e. g. in hospitals, old people's homes, nursing homes etc. and doctor's surgeries for the collection of body fluids or to drop instruments according to their use and before their cleaning.

Restriction of use:

Hot and/or frozen substances, alcoholic solutions, microwave, sterilization processes. Do not reuse!

Storage and transport conditions:

Product should be kept in a dry environment. If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).

Other:

Goods are sensitive to mechanical impact and must not be subjected to impacts or high pressure be exposed. Goods must be protected from moisture and humidity. Use outside the intended purpose is the responsibility of the user.

Version 13.06.2023



Omni-Pac Group
Second Avenue,
Flixborough Industrial Estate
Scunthorpe
North Lincolnshire DN15 8SD
UNITED KINGDOM
Tél: +44 (0)1724 847 207
www.omnipacgroup.co.uk

