

Medical CarePulp® : Versatile Solutions

GP BOWL



1L

Registered medical product with 

Dimensions piece 195mm x 195mm x 82mm (L x W x D)

Weight / 22 g

piece Carton

Weight of the Carton : 4.4 kg

Carton dimensions: 195mm x 195mm x 82mm (L x W x D)

Number of pieces per carton : 200

Number of carton per palette : 24

Number of pieces per pallet : 4,800

Pallet dimensions: 1200mm 1000mm x 1370mm (L x W x H)

No printing

Available colors

Grey

UDI-DI 5060802530032

Product should be kept in a dry environment.
If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).
The characteristics are given as an indication. Non contractual document. Changes excepted. Weights & measures are of average data!
You can find more information about this product on the following page

Version 13.06.2023



Omni-Pac Group
Second Avenue,
Flixborough Industrial Estate
Scunthorpe
North Lincolnshire DN15 8SD
UNITED KINGDOM
Tél: +44 (0)1724 847 207
www.omnipacgroup.co.uk



Supplement to the product data sheet

Medical products

Product description

Composition (materials used):

Moulded fibre products made of recycled paper and fresh fibre or only recycled paper with the addition of process-related and functional additives.

Intended use:

As medical devices for use in the medical field, e. g. in hospitals, old people's homes, nursing homes etc. and doctor's surgeries for the collection of body fluids or to drop instruments according to their use and before their cleaning.

Restriction of use:

Hot and/or frozen substances, alcoholic solutions, microwave, sterilization processes. Do not reuse!

Storage and transport conditions:

Product should be kept in a dry environment. If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).

Other:

Goods are sensitive to mechanical impact and must not be subjected to impacts or high pressure be exposed. Goods must be protected from moisture and humidity. Use outside the intended purpose is the responsibility of the user.

Version 13.06.2023



Omni-Pac Group
Second Avenue,
Flixborough Industrial Estate
Scunthorpe
North Lincolnshire DN15 8SD
UNITED KINGDOM
Tél: +44 (0)1724 847 207
www.omnipacgroup.co.uk

