Medical CarePulp[®]: Versatile Solutions Male Urinal Bottle Round

Registered medical product with $\mathbf{C}\mathbf{\epsilon}$

OMNI-PAC



875 ML

Dimensions piece	260 x 100 x 130 mm (L x W x D)	
Weight /	30 g	
piece Carton	Dimensions: 260 x 100 x 130 mm (L x W x D) Weight per piece: 30 g Maximum volume 875 ml Number of pieces per carton: 100 Storage: Cool, dry conditions Shelf life: 3 years	
No printing		
Available colors		
Grey	UDI-DI 5060802530001	
The characteristics are given as You can find more information a	an indication. Non contractual document. Changes excepted. Weights & measures are of average data! bout this product on the following page	Version 13.06.2023
		-
Č	Omni-Pac Group Second Avenue, Flixborough Industrial Estate Scunthorpe North Lincolnshire DN15 8SD	OMNI-PAC

Tél: +44 (0)1724 847 207 www.omnipacgroup.co.uk

UNITED KINGDOM

ACT NOW!



Supplement to the product data sheet



Medical products

Product description

Composition (materials used):

Moulded fibre products made of recycled paper and fresh fibre or only recycled paper with the addition of process-related and functional additives.

Intended use:

As medical devices for use in the medical field, e. g. in hospitals, old people's homes, nursing homes etc. and doctor's surgeries for the collection of body fluids or to drop instruments according to their use and before their cleaning.

Restriction of use:

Hot and/or frozen substances, alcoholic solutions, microwave, sterilization processes. Do not reuse!

Storage and transport conditions:

Product should be kept in a dry environment. If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).

Other:

Goods are sensitive to mechanical impact and must not be subjected to impacts or high pressure be exposed. Goods must be protected from moisture and humidity. Use outside the intended pur-pose is the responsibility of the user.



Omni-Pac Group Second Avenue, Flixborough Industrial Estate Scunthorpe North Lincolnshire DN15 8SD UNITED KINGDOM Tél: +44 (0)1724 847 207 www.omnipacgroup.co.uk

