# Medical $CarePulp^{\text{\tiny B}} \ : \text{Versatile Solutions}$ Midi Slipper Pan

Registered medical product with  $\mathbf{C}\mathbf{\epsilon}$ 



1.3LTR

# **Dimensions piece** 40 x 300 x 70 mm (L x W x D) Weight 44 g 1 Dimensions: 340 x 300 x 70 mm (L x W x D) piece Carton Weight per piece: 44 g Maximum volume 1.5 L Number of pieces per carton: 200 Storage: Cool, dry conditions Shelf life: 3 years No printing Available colors Grey UDI-DI 060802530070

The characteristics are given as an indication. Non contractual document. Changes excepted. Weights & measures are of average data! You can find more information about this product on the following page



**Omni-Pac Group** Second Avenue, Flixborough Industrial Estate Scunthorpe North Lincolnshire DN15 8SD UNITED KINGDOM Tél: +44 (0)1724 847 207 www.omnipacgroup.co.uk



Version 13.06.2023

## Supplement to the product data sheet



## Medical products

### **Product description**

Composition (materials used):

Moulded fibre products made of recycled paper and fresh fibre or only recycled paper with the addition of process-related and functional additives.

### Intended use:

As medical devices for use in the medical field, e. g. in hospitals, old people's homes, nursing homes etc. and doctor's surgeries for the collection of body fluids or to drop instruments according to their use and before their cleaning.

### Restriction of use:

Hot and/or frozen substances, alcoholic solutions, microwave, sterilization processes. Do not reuse!

Storage and transport conditions:

Product should be kept in a dry environment. If stored in a non airconditioned environment, please protect the products from direct humidity and from too humid and too dry air (blowpipe).

#### Other:

Goods are sensitive to mechanical impact and must not be subjected to impacts or high pressure be exposed. Goods must be protected from moisture and humidity. Use outside the intended pur-pose is the responsibility of the user.



Omni-Pac Group Second Avenue, Flixborough Industrial Estate Scunthorpe North Lincolnshire DN15 8SD UNITED KINGDOM Tél: +44 (0)1724 847 207 www.omnipacgroup.co.uk

